

Used 15,000-Gallon Porta-Pac

Year	Capacity	PSI	Serial No. / National Board No.	U1A Data Report
1969	15,000	250	370576 / N.A.	N.A.

Porta-Pac is equipped with (4) Emergency Shutoff Boxes and ready for nitrogen bottles for ESS hookup. Includes (2) 250PSI Reliefs, Visible Liquid Level, Temperature and Pressure Gauges.

Liquid In: 2" MEC Internal Valve w/ (2) 2" Sight Flow Swing Checks, (2) 2" Angle Valves, (2) 2" Break-away and 3.25" acme for transport hook up.

Vapor Equalization: 1.25" MEC Internal valve w/ (2) 1.25" MEC ESV's, (2) 1.25" Angle valves, (2) 1.25" Break-aways and 1.75" acme for transport hook up.

Liquid Out: 3" MEC Internal Safety valve with 2" 300# MEC Globe valve. Customer to make final connection for service.

